

notice the angle

# **HyProCure<sup>TM</sup>** Sinus Tarsi Implant

not “Just Another Subtalar Arthroereisis Implant”

**When orthotics are ineffective**

**or**

**Rearfoot reconstructive surgery  
is too aggressive**

**Think HyProCure<sup>TM</sup>**

\* **Fills in the sinus and canalis portion of the sinus tarsi**


\* **Self-centering / Self-seating design**

\* **No bone cutting or drilling**

\* **Commonly used in pediatric, adult,  
and geriatric patients**

\* **Completely reversible**

\* **Minimally invasive**



**Before**



**After**



**Before**



**After**

# Surgical Technique



An oblique linear skin incision is made within the skin tension lines 1.5 to 2.0 cm long over the sinus tarsi.



The incision is bluntly dissected with curved hemostats through the subcutaneous tissue to the anteriolateral aspect of the sinus tarsi using tenotomy scissors, the interosseous ligament is transected. A guide wire is placed into the sinus tarsi.



The number six trial sizer is inserted into the sinus tarsi from anteriolateral to posteromedial and the foot is maximally pronated. The goal is to restore the normal amount of pronation. Use larger trial sizers sequentially, until the desired size is determined.



The sizer is removed and the appropriate stent is selected, placed on the driver, and advanced into the sinus tarsi.



The driver is then removed and the incision is irrigated. The incision is closed per the surgeon's choice and a light compressive dressing is applied.

## Indications

- \* **Obliteration (partial/full) of the sinus tarsi**
- \* **Flexible and reducible deformity**
- \* **Talar deviation medial, and/or anterior, and/or plantar on the calcaneus**

## Contra-Indications

- \* **Rigid deformity - non reducible**
- \* **Extremely poor bone stock**
- \* **Other conditions that may place the patient at risk**

Please visit **[www.HyProCure.com](http://www.HyProCure.com)** for videos and more information, including:

what is hyperpronation?

common symptoms

why HyProCure™?

testimonials

physician information

frequently asked questions

news & events



Distribütör :



[www.integramedikal.com](http://www.integramedikal.com)  
[info@integramedikal.com](mailto:info@integramedikal.com)

Tel : 0(216) 326 67 25 Fax : 0(216) 546 07 08

Bölge Bayisi :



ISO 13485:2003  
 © 2008 GraMedica  
 Rev 11.12.08